

March 22, 2004


 HealWell
Buy  Online!
[Home](#) | [Contact Us](#) | [Email Us](#)
[About Us](#) | [Find a Retailer Near You](#)

The liners of the Clinitek® HealWell™ Multi AFO/Contracture Splints with Microban® are treated with Microban® Antimicrobial Protection.

WHAT IS



Microban® is very effective at inhibiting microbial growth. Compare growth around Microban® treated orthopedic fabric (right) to non-treated sample (left). A clean zone of inhibition is observed with the Microban® treated sample.

Microban® is an antimicrobial agent which, when introduced into the actual molecular structure of fibers, fabrics and polymers, provides continuous inherent control over the growth of certain odor and stain causing bacteria, yeasts and fungi.

Contracture splints can become a hostile healing environment. The combination of warmth and moisture in the contracture splint liner forms an ideal breeding ground for certain bacteria, yeasts and fungi. **Microban® cannot prevent infection, but it can help maintain a fresher product in the context of a good patient hygiene program.**



The Microban® protected liners can prevent the growth of certain odor and stain causing bacteria, yeast and fungi. The result is a more hygienic, fresher product for improved patient compliance and comfort. **Microban® is not intended to prevent infection.**

Toxicity and dermal studies have shown products with Microban® to be completely safe, non-toxic, and non-irritating to the skin. Microban® is effectively and safely used currently in a number of medical and consumer products.

Copyright © 2001 - 2004 [FLA Orthopedics, Inc.](#)
All Rights Reserved

Visit some of our other sites.



